

Gas drying and conditioning system

EmiTest GUARD

Description

EmiTest *Guard* is made to remove easily and efficiently water in order to protect samplers such as EmiTest *Iso*.


Its main features are:

- Easy to handle – all-in-one solution
- All inclusive system between probe and sampler
- Easy tightening
- Minimises the use of glassware
- Water condensate collection in standard laboratory bottle
- Integrated silica gel trap



Principle

- ① Sampled gas is cooled in a serpentine contained in the cooling bath (possibility to use ice or cold water to increase efficiency)
- ② Condensates are collected in a dedicated bottle
- ③ Remaining humidity is removed by silica gel trap



Technical characteristics

General specifications	
Dimensions (LxWxH)	340x222x570 mm
Weight	8 kg
Reference	ES300-0000

Technical specifications may change without previous warning